



DIAGNOSTIC  
LABORATORY  
OF OKLAHOMA.

# Simplifying STAT Testing

DLO is simplifying our STAT process with a standardized testing menu and a single drop off location

## DLO STAT Testing Offering for OKC

DLO is pleased to offer a simplified and standard STAT testing menu available at the DLO Central Lab, located at 225 NE 97th Street in Oklahoma City. DLO will no longer accept STAT drop offs at the INTEGRIS hospital labs in the Oklahoma City metro.

STAT testing results will be available within two hours of the specimen drop off time at this location. Please see the reverse side for DLO's complete standard STAT menu list and specimen requirements.

If you require specimen transportation, DLO Logistics will be available for STAT pick-ups within the STAT service area (see map). STAT test results will be available four hours after the pick-up request has been received when using DLO Logistics. To request a STAT pick-up, call 405.608.6100, option 3, or place your request online through Quamum HCP at [dlolab.com/pickup](http://dlolab.com/pickup).


If multiple tests are ordered on a single requisition, you must specify which test(s) need to be run STAT.

As a reminder, most in-house non-STAT testing is resulted the following day by 8 a.m., and STAT specimens may not be dropped off at DLO Patient Service Centers.

If you require tests outside of the standard STAT menu on a consistent basis, contact your DLO Account Executive to see if a STAT Service Agreement can be completed.

Diagnostic Laboratory of Oklahoma, DLO, any associated logos, and all associated Diagnostic Laboratory of Oklahoma registered or unregistered trademarks are the property of Diagnostic Laboratory of Oklahoma. © 2023 Diagnostic Laboratory of Oklahoma, L.L.C. All rights reserved. 6/2023

### DLO Oklahoma City STAT Logistics Service Area



DIAGNOSTIC  
LABORATORY  
OF OKLAHOMA.

800.891.2917 • [dlolab.com](http://dlolab.com)

# Tests Eligible for STAT Testing and Specimen Collection Requirements

TEST CODE	TEST NAME	DLO CENTRAL LAB SPECIMEN TYPE
823	ALANINE TRANSAMINASE	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
223	ALBUMIN	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
234	ALKALINE PHOSPHATASE	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
5509	AMMONIA (P)	EDTA PLASMA, FROZEN
243	AMYLASE	SERUM, SST
822	ASPARTATE TRANSAMINASE	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
10165	BASIC METABOLIC PANEL	SERUM SST
285, 34162	BILIRUBIN, DIRECT	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER- CENTRIFUGE/LIGHT PROTECT
6631	BILIRUBIN, FRACTIONATED,	AMBER MICROTAINER SST X 2 OR LITHIUM HEPARIN WITH GEL
287	BILIRUBIN, TOTAL	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER- CENTRIFUGE/LIGHT PROTECT
20381		
37386 or 59630		
303		
310		
6399		
20253	BILIRUBIN, TOTAL,	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
330	B-TYPE NATRIURETIC	EDTA PLASMA, FROZEN
396		
10231	CALCIUM	SERUM OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
97637	CARBON DIOXIDE	SERUM OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
97638		
97636	CBC (INCLUDES DIFF/PLT)	EDTA WHOLE BLOOD
374		
375	CBC (INCLUDES DIFF/PLT)	EDTA WHOLE BLOOD
4420		
418	CHLORIDE	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
34392	CHORIONIC GONADOTROPIN,	FRESH, CLEAN CATCH IN STERILE CONTAINER OR STOCKWELL TUBE
443		
483	COMPREHENSIVE METABOLIC	SERUM SST
8435	COVID, FLU A/B PANEL***	ANTERIOR NARES/NASAL (AN) OR NASOPHARYNGEAL (NP) PHOSPHATE BUFFED SALINE (PBS), VCM, UTM
8396	COVID, RSV PANEL ***	NASOPHARYNGEAL (NP) PHOSPHATE BUFFED SALINE (PBS), VCM, UTM
509		
510	COVID, FLU A/B, AND RSV***	NASOPHARYNGEAL (NP) PHOSPHATE BUFFED SALINE (PBS), VCM, UTM
7998		
1759	CREATINE KINASE, TOTAL	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
10256	CREATININE W/EGFR	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
654	CRP	SERUM WITH OR WITHOUT GEL BARRIER
11177		
593	DIGOXIN	SERUM WITHOUT GEL BARRIER
606	ELECTROLYTE PANEL	SERUM SST
622	ETOH	SERUM WITH OR WITHOUT GEL BARRIER
763		
718	GLUCOSE, SERUM	SPUN SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
733	HCG, TOTAL, QL	SERUM SST
754	HCG, TOTAL, QN	SERUM SST
8847		
10314	HEMATOCRIT	EDTA WHOLE BLOOD
10553	HEMOGLOBIN	EDTA WHOLE BLOOD
59039	HEMOGLOBIN & HEMATOCRIT	EDTA WHOLE BLOOD
905		
3020	HEMOGRAM/PLT	EDTA WHOLE BLOOD
5463	HEPATIC FUNCTION PANEL	SERUM SST
7909	HETEROPHILE, MONO SCREEN	SERUM SST
	INFLUENZA A AND B ANTIGEN	N/P SWAB IN VCM
	LDH, TOTAL	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
	LIPASE	SERUM SST
	MAGNESIUM	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
	PARTIAL THROMBOPLASTIN	SODIUM CITRATE
	PHOSPHORUS, INORGANIC	SERUM SST
	POTASSIUM	SERUM SST
	PROTEIN, TOTAL	SERUM SST
	PROTHROMBIN TIME WITH INR	SODIUM CITRATE
	RENAL FUNCTION PANEL	SERUM SST
916	STREPTOCOCCUS GROUP A AG	DUAL SWAB RED CAP AMIES
8578	TROPONIN I	SERUM OR LITHIUM HEPARIN PLASMA, WITH OR WITHOUT GEL BARRIER
917	URIC ACID	SERUM SST OR LITHIUM HEPARIN PLASMA, WITH GEL BARRIER
8579	URINALYSIS, COMPLETE	FRESH, CLEAN CATCH IN STOCKWELL TUBE AND URINE CULTURE TUBE (GRAY-TOP)
	URINALYSIS, COMPLETE	FRESH, CLEAN CATCH IN STOCKWELL TUBE
	URINALYSIS, REFLEX	FRESH, CLEAN CATCH IN STOCKWELL TUBE
	VALPROIC ACID	SERUM WITHOUT GEL BARRIER

\*\*\* ONLY the FLU A/B and RSV testing will be performed STAT. COVID testing is not performed STAT. Turn-around-time for COVID testing is 24 hours from time of receipt.

